

## Produktinformation DHS- Platten / Product information DHS Hip screw Plates



Material / Material :	ISO 5832-1
Schaftdicke / Thickness	6,0mm
Schaftbreite / Width :	19,0mm
Lochabstand / Hole spacing:	16,0mm
Zylinderlänge / Barrel length:	38,0mm

DHS-Hüftplatte 135° 2-Loch, Länge = 46mm, unsteril DHS Hip screw Plate 135° 2-Hole, Length = 46mm, non-sterile
DHS-Hüftplatte 135° 3-Loch, Länge = 62mm, unsteril DHS Hip screw Plate 135° 3-Hole, Length = 62mm, non-sterile
DHS-Hüftplatte 135° 4-Loch, Länge = 78mm, unsteril DHS Hip screw Plate 135° 4-Hole, Length = 78mm, non-sterile
DHS-Hüftplatte 135° 5-Loch, Länge = 94mm, unsteril DHS Hip screw Plate 135° 5-Hole, Length = 94mm, non-sterile
DHS-Hüftplatte 135° 6-Loch, Länge = 110mm, unsteril DHS Hip screw Plate 135° 6-Hole, Length = 110mm, non-sterile
DHS-Hüftplatte 135° 7-Loch, Länge = 126mm, unsteril DHS Hip screw Plate 135° 7-Hole, Length = 126mm, non-sterile
DHS-Hüftplatte 135° 8-Loch, Länge = 142mm, unsteril DHS Hip screw Plate 135° 8-Loch, Length = 142mm, non-sterile
DHS-Hüftplatte 135° 10-Loch, Länge = 174mm, unsteril DHS Hip screw Plate 135° 10-Hole, Length = 174mm, non-sterile
DHS-Hüftplatte 135° 12-Loch, Länge = 206mm, unsteril DHS Hip screw Plate 135° 12-Hole, Length = 206mm, non-sterile
DHS-Hüftplatte 135° 14-Loch, Länge = 238mm, unsteril DHS Hip screw Plate 135° 14-Hole, Length = 238mm, non-sterile
DHS-Hüftplatte 135° 16-Loch, Länge = 270mm, unsteril DHS Hip screw Plate 135° 16-Hole, Length = 270mm, non-sterile