

ZWEYMÜLLER*Type

Femoral Stem

Based on the design principles of the well known and clinically proven press-fit Femoral Stem, conceived by Professor Karl Zweymüller in Austria in 1979

* Zweymüller, K. (1986)
A cemented titanium hip endoprosthesis system, based on press-fit fixation: Basic research and clinical results (Inst. Course lect., 35:203-205)

Available in 11 sizes as
Titanium Forgings or Finished Stems
Standard or Lateral

FORGED IN
GERMANY



Your Contract Manufacturer
for Generic Hip Systems

Aristotech[®]
Orthopaedic Forgings and Services

ZWEYMÜLLER

Type

Femoral Stem

Titanium

Forgings



11

10

09

08

07

06

05

04

03

02

01

standard & lateral

Aristotech®

2 CCD
angles
123°/131°

Finished Stems



01

02

03

04

05

06

07

08

09

10

11

standard & lateral

FRÜHGE

Datum

ARISTOTECH®
Digital Technologies

Zul. Abw. / Oberfläche

Entform
Verset
Durchbieg
Verwerfu
Maße
werk
Flach

Aristotech® Orthopaedic Forgings and Services




Zweymüller Type

Femoral Stems and Rasps

Order Information - Catalog Numbers

	Forgings Titanium	Finished Stems	Finished Stems		
size	blanc	standard	lateral	size	Rasp
01	A1201.01-03	A1201.01	A1202.01	01	AT71201
02	A1201.02-03	A1201.02	A1202.02	02	AT71202
03	A1201.03-03	A1201.03	A1202.03	03	AT71203
04	A1201.04-03	A1201.04	A1202.04	04	AT71204
05	A1201.05-03	A1201.05	A1202.05	05	AT71205
06	A1201.06-03	A1201.06	A1202.06	06	AT71206
07	A1201.07-03	A1201.07	A1202.07	07	AT71207
08	A1201.08-03	A1201.08	A1202.08	08	AT71208
09	A1201.09-03	A1201.09	A1202.09	09	AT71209
10	A1201.10-03	A1201.10	A1202.10	10	AT71210
11	A1201.11-03	A1201.11	A1202.11	11	AT71211



Handle AT176000



additional
Generic Femoral Hip Systems
available

Aristotech®
Your Orthopaedic
Contract Manufacturer



Wagner Revision
Type



FemoraLock
Style



Spotorno
Type



Müller
Type



BioRail
Style



Dual Mobility Cup



PoroCup

Forgings or Finished
at **NO TOOLING COSTS**

Ask for more Information

ARISTOTECH®
Implant Technologies

Zul. Abw. / Oberfläche

Maßstäbe

Werkstoff

Wahlver

Ultra Modern Facility in Luckenwalde near Berlin / Germany



Founded in 2007, Aristotech manufactures for OEM Orthopaedic Companies all over the world cost-effective Hip Replacement Forgings and finished Implants in a broad range of capabilities with validated processes and services.

- Net Shape Forgings - Precision Machining - Cleaning and Passivation -
- Clean Room Assembly - Sterile Packing Services -



Supporting especially small and medium sized
OEM Orthopaedic Companies
with Net Shape Forgings or Finished
Cost-Effective Generic Hip Systems
At No Tooling Costs

www.aristotech.de

D-14943 Luckenwalde +49 3371 40640-0

service@aristotech.de