

WAGNER* Type

Femoral Stem

Taped fluted Titanium Femoral Stem

for revision or complex primary THA
based on the design of
Professor Heinz Wagner
developed in Germany in 1987

*Reference : Wagner H.


A revision prosthesis
for the hip joint.

Orthopaede 1989, 18. 438-453

Available as:

Titanium Forgings or finished Stems
in 4 lengths and 12 diameters

FORGED IN
GERMANY



Your Contract Manufacturer
for Generic Hip Systems

Aristotech
Orthopaedic Forgings and Services



WAGNER Type

Femoral
Stem

305 mm

265 mm

225 mm

190 mm

finished

Aristotech®

WAGNER Type finished

Order Information - Catalog Numbers

Reamer	diameters in mm	Stem length 190 mm	Stem length 225 mm	Stem length 265 mm	Stem length 305 mm
A4094.80	Ø 14	A4087.14	A4086.14	A4085.14	A4084.14
A4094.81	Ø 15	A4087.15	A4086.15	A4085.15	A4084.15
A4094.82	Ø 16	A4087.16	A4086.16	A4085.16	A4084.16
A4094.83	Ø 17	A4087.17	A4086.17	A4085.17	A4084.17
A4094.84	Ø 18	A4087.18	A4086.18	A4085.18	A4084.18
A4094.85	Ø 19	A4087.19	A4086.19	A4085.19	A4084.19
A4094.86	Ø 20	A4087.20	A4086.20	A4085.20	A4084.20
A4094.87	Ø 21		A4086.21	A4085.21	A4084.21
A4094.88	Ø 22		A4086.22	A4085.22	A4084.22
A4094.89	Ø 23			A4085.23	A4084.23
A4094.90	Ø 24			A4085.24	A4084.24
A4094.91	Ø 25			A4085.25	A4084.25





Ø 25 mm



Ø 21 mm



Ø 17 mm

WAGNER Type

Blanc
Forgings
Titanium

forged

Aristotech®

WAGNER Type Blanc Forgings Titanium

Order Information - Catalog Numbers

diameters in mm	Forging blanc 190 mm	Forging blanc 225 mm	Forging blanc 265 mm	Forging blanc 305 mm
Ø 14	A4094.50-03	A4094.50-03	A4094.50-03	A4094.50-03
Ø 15	A4094.50-03	A4094.50-03	A4094.50-03	A4094.50-03
Ø 16	A4094.50-03	A4094.50-03	A4094.50-03	A4094.50-03
Ø 17	A4094.50-03	A4094.50-03	A4094.50-03	A4094.50-03
Ø 18	A4094.60-03	A4094.60-03	A4094.60-03	A4094.60-03
Ø 19	A4094.60-03	A4094.60-03	A4094.60-03	A4094.60-03
Ø 20	A4094.60-03	A4094.60-03	A4094.60-03	A4094.60-03
Ø 21	A4094.60-03	A4094.60-03	A4094.60-03	A4094.60-03
Ø 22	A4094.70-03	A4094.70-03	A4094.70-03	A4094.70-03
Ø 23	A4094.70-03	A4094.70-03	A4094.70-03	A4094.70-03
Ø 24	A4094.70-03	A4094.70-03	A4094.70-03	A4094.70-03
Ø 25	A4094.70-03	A4094.70-03	A4094.70-03	A4094.70-03



Schnitt A-A

Schnitt C-C

PRÜFGABE

Datum

ARISTOTECH®
Innov. Technologie

Zul. Abw. / Oberfläche

Höchst

Werkst

Richtst

additional
Generic Femoral Hip Systems
available


Aristotech[®]
Your Orthopaedic
Contract Manufacturer



BioRail
Style



FemoraLock
Style



Spotorno
Type



Müller
Type



Zweymüller
Type


Dual Mobility Cup



PoroCup

Forgings or Finished
at **NO TOOLING COSTS**

Ask for more Information



Ultra Modern Facility in Luckenwalde near Berlin / Germany



Founded in 2007, Aristotech manufactures for OEM Orthopaedic Companies all over the world cost-effective Hip Replacement Forgings and finished Implants in a broad range of capabilities with validated processes and services.

- Net Shape Forgings - Precision Machining - Cleaning and Passivation -
- Clean Room Assembly - Sterile Packing Services -



Supporting especially small and medium sized
OEM Orthopaedic Companies
with Net Shape Forgings or Finished
Cost-Effective Generic Hip Systems
At No Tooling Costs

www.aristotech.de

D-14943 Luckenwalde +49 3371 40640-0

service@aristotech.de